

Low-cost, portable, easy-to-use kiosks to facilitate home-cage testing of non-human primates during vision-based behavioral tasks

Christopher Giverin, Hamidreza Ramezanpour, Kohitij Kar

Department of Biology, Centre for Vision Research, York University, Toronto, ON, Canada
Contact: koh1t1j@yorku.ca, cgiverin@yorku.ca

- 1 Non-human primates, especially rhesus macaques, have played a significant role in our current understanding of the neural computations
2 underlying human vision. Apart from the established homologies in the visual brain areas between these two species, and our extended
3 abilities to probe detailed neural mechanisms in monkeys at multiple scales, one major factor that makes NHPs an extremely appealing animal
4 model of human-vision is their ability to perform human-like visual behavior. Traditionally such behavioral studies have been conducted in
5 the controlled laboratory settings. Such in-lab studies offer the experimenter a tight control over many experimental variables. However,
6 typical constraints related to such experiments include, 1) limited total experimental time, 2) requirement of dedicated human experimenters
7 for the NHPs, 3) requirement of additional lab-space, 4) additional time and training required for chairing and head restraints of monkeys,
8 5) NHPs often need to undergo invasive surgeries for a head-post implant. To overcome these limitations, many laboratories have now
9 adapted home-cage behavioral training and testing of NHPs. Home-cage behavioral testing enables the administering of many vision-based
10 behavioral tasks simultaneously across multiple monkeys with much reduced human personnel requirements, no NHP head restraint, and
11 provide NHPs access to the experiments without specific time constraints. Here, we provide the details of operating and building a portable,
12 easy-to-use kiosk for conducting home-cage vision-based behavioral tasks in NHPs.

Non-human primates | kiosk | home-cage | visual psychophysics

1 SUMMARY

- 2 This article will summarize the operation and building blocks of a
3 low-cost, portable, easy-to-use kiosk developed in our laboratory
4 at York University (Toronto, Canada). The design and motivation

5 for the kiosk have been heavily inspired by earlier work from Elias
6 Issa in the laboratory of James J DiCarlo at the McGovern Institute
7 of Brain Research at MIT ([Rajalingham et al., 2018](#); [Kar et al., 2019](#)).
8 The primary goal of this article is to make the details of our
9 version of the kiosks readily available across the NHP laboratories
10 worldwide. An elaborate pictorial depiction of the kiosk (as provided
11 here) will also allow laboratories (that have mechanical workshops




Fig. 1. Lifting the kiosk. Primarily, the kiosk is lifted by the experimenter immediately before placing it on the cage. There are two handles of the side that allow the experimenter to lift the ~25-30 lbs kiosk.

12 available to them via their university or departmental resources)
 13 to get a gist of the mechanisms involved and quickly produce a
 14 customized version for use. Adopting these kiosks will enhance
 15 NHP research, and openly sharing our design will help us improve
 16 them in the future. We also expect to update this article as we
 17 continue incorporating additional capabilities into this system.

18 Interaction of the human experimenter with the kiosk

19 The kiosk is approximately 25-30 lbs. in weight. The National
 20 Institute for Occupational Safety and Health (NIOSH) has created
 21 a mathematical model that can predict the chance of injury
 22 when lifting weights, taking into consideration factors such as
 23 weight, frequency of lifting, and movement during lifting. A lifting

24 equation (see <http://www.cdc.gov/niosh/docs/94-110/>), calculates
 25 a recommended weight limit for one individual in various situations.
 26 The lifting equation establishes a maximum load of 51 pounds,
 27 which is then adjusted based on various factors, including how
 28 often the weight is lifted, the twisting of the back during lifting, the
 29 vertical distance the load is lifted, the distance of the load from the
 30 body, the distance moved while lifting the load, and the ease of
 31 holding onto the load. Given that the kiosk can be carried around
 32 in a cart, and only needs to be lifted right before placing it on the
 33 cage bars, the overall time of lifting as well as its total weight fits
 34 the capacity of most experimenters and are within the NIOSH
 35 estimated limits. **The experimenter should lift the kiosk by its**
 36 **side using the protruding handles (Figure 1).** In an alternate
 37 model of the kiosk, the handles can be also placed vertically
 38 (instead of horizontally).
 39
 40
 41




Fig. 2. Interacting with the kiosk. **A.** The top panel that remains attached to the kiosk can be opened to access the rest of the kiosk. **B.** One of the detachable inserter panels holds the tablet/touch screen. **C.** A protective front panel is inserted at the end of the experimental session to protect the electronics at the back of the kiosk from the monkey. It is also removed at the beginning of the session to allow the monkey full access to the touch screen. **D.** A removable bottle is placed in a dedicated attachment in the back of the kiosk to hold the reward fluid (typically water). The bottle can be removed daily to refill the fluid. **E.** The back door (that stays attached to the kiosk) can be opened to insert the electronics, and make proper adjustments (e.g., tighten the tablet into the inserter).




Fig. 3. The touch screen/tablet inserter panel. **A.** *Left Panel:* The view from the experimenter's end. The knobs can be loosened to insert the tablet from the top of the kiosk. The back of an iPad is visible in the figure as an example configuration. *Right Panel:* The view from the monkey's end. The visual stimulus can be seen in the example. The monkeys can touch this screen to indicate their choices. **B.** Possible items to add to the setup via this inserter. Apertures can be drilled at the bottom to allow for two IR illuminators. A small hole at the top of the inserter panel can be used to insert a high-definition USB camera with an IR filter. These can be combined to develop a low-cost eye tracking system (available in future versions of the kiosk).




Fig. 4. The dimensions of the touch screen/tablet inserter panel. The exact dimensions can vary based on the make of the tablet used. The aperture size can be smaller than the touch screen if the entire screen is not utilized in the tasks. However, the size shouldn't be larger than the touch screen since that will allow the monkeys to access the electronics (usually placed behind the touch screen).




Fig. 5. Securely attaching the kiosk to the cage. **A.** The kiosk can be attached to the cage **B.** The top part of the kiosk slides into the cage, while the bottom has a locking mechanism that holds it in place such that the NHP cannot move the kiosk. **C.** Once placed correctly has no exposed gaps that can encourage the NHP to shake the kiosk. To get a closer look at how the kiosk should be placed, kindly refer to the following video, <https://github.com/vital-kolab/nhp-turk/blob/main/videos/Placement.mp4>

42 Attaching the kiosk to the cage

43 The kiosk needs to be securely attached to the top (**Figure 5B:**
44 top panel) and the bottom (**Figure 5B:** bottom) of the cage opening
45 with specific mechanisms (see our demo videos on GitHub
46 [Placement.mp4](#) for details and part list). Please refer to the sup-
47 plemental video to observe the exact locking mechanisms that
48 prevent the kiosk from being detached from the cage (regardless
49 of any force applied by the NHP). Given that NHP cages come in
50 various configurations, the exact attachment styles will likely need
51 to be modified to cater to individual caging systems.

52 The internal components of the kiosk

53 When the monkey has access to the kiosk (that is when the kiosk
54 is securely attached to the cage), experimenters can still access
55 parts of the kiosk to make necessary changes to the experiments,
56 change out the reward fluid, replace batteries for the tablet (or other
57 electronics, e.g., Arduinos, Raspberry Pi). **Majority of the com-**
58 ponents like the side walls of the kiosk, the inserter panels,
59 and the top and the back panel doors are all polycarbonate
60 sheets. All parts have been itemized (see e.g. in **Figure 6**) and
61 shared at [Bill of materials.xlsx](#). The following are the most relevant
62 components of the kiosks with respect to the experimenter's needs:

63 **The top panel (non-detachable from the kiosk).** The top panel
64 door can be pulled open to access the content of the kiosk (**Fig.**
65 **2A**). Typically, one would open it in the beginning of the experimen-
66 tal session to insert the tablet/touch screen into the kiosk. The
67 experimenter can also place dry treats and toys into the front end
68 of the kiosk during initial days of the NHP training to make the NHP
69 feel motivated to insert their hands into the kiosk. They can also
70 adjust the electronics at the back of the kiosk.

71 **The two inserter panels (removable from the kiosk).** The kiosk
72 uses two panels that can be inserted into it. The front panel (**Fig.**

73 **2C**) is for protective purposes, and it is only inserted at the end of
74 the experimental session, so that the monkey cannot access the
75 electronics at the back of the kiosk.

Example contents of the Bill of materials


Part No. 1

Part Name	Polycarbonate Sheet (Opaque)
Description	Inside dividers and Tablet inserters
Vendor	McMaster - Carr
Vendor Part No.	8574K211
Quantity	8
Size per unit	12" x 12" x 1/4"
Unit cost (USD)	16.94
Total cost (USD)	135.52
URL	https://www.mcmaster.com/8574K211/


Fig. 6. Example contents of the Bill of materials. Please find an excel sheet with detailed information on all of these items for each part at <https://github.com/vital-kolab/nhp-turk/blob/main/Bill.of.Materials.xlsx>

76 Typically the tablet is removed for charging at the end of the
77 day. Therefore the front panel has to be placed. This panel is also

A Kiosk Cage assembly **B** Outside the cage **Inside the cage**



C Front of the kiosk facing the NHP



D Sideview of the kiosk




Fig. 7. Kiosk components critical for NHP interactions. **A.** Kiosk Cage assembly. Before the NHP is allowed to interact with the tablet by removing the front-most protective inserter panel, the kiosk has to be securely attached to the cage **B.** The monkey accesses the touch screen by sitting right in front of the kiosk on a horizontal cage divider panel. **C.** The front of the kiosk facing the NHP has two apertures. The top one allows for clear visual access to the NHP. The bottom one allows the NHP to insert its hand and interact with the touch screen. **D.** Sideview of the kiosk illustrates the space at the front that can be used to insert toys and treats to initially encourage the NHP to get their hands in via the aperture.

removed at the beginning of the experimental session so that the NHP can access the touch screen. The second inserter panel (**Fig. 2B**) holds the tablet/touch screen. This panel, ideally could be left in the kiosk at all times. The back of this panel (that is accessible to the experimenter when the NHP has access to the kiosks), has brackets and knobs to securely place an iPad or other tablets into that slot (see **Figure 3A; left panel**). The front of this panel has an aperture sized according to the dimension of the touch screen (**Figure 3A; right panel**). Depending on whether the vision tasks are conducted on iPads or Android tablets or via other means, the size of the aperture should be adjusted. The current dimensions of this inserter panel is provided in **Figure 4**. Other features to the kiosks, like eye tracking can be added by making modifications to this inserter panel (**Figure 3B**). For instance, an example strategy could involve making three additional holes on this inserter panel. Two of these could be at the bottom to accommodate two IR illuminators and one could be at the top to allow for a HD USB-

based camera with IR filters. A combination of these can be used to set up a low cost eye tracker (*more details of this tracker will be published soon*).

The back panel (non-detachable from the kiosk). The back panel allows the experimenter to insert the electronics (**Fig. 2E**) associated with running the experiments. It is also attached to the section that contains the reward fluid bottle (**Fig. 2D**). The bottle has a tube inserted in it that is usually connected to a pump and then via another tube to a front nozzle (that the NHP uses to get the reward).

Interaction of the animal with the kiosk

To begin the experimental session, the kiosk must be first securely attached to the cage (**Figure 7A**). The kiosk can be placed both at the top-half or the bottom half of a standard NHP caging system. If it is placed in the top-half, a horizontal divider should be in place

110 to allow the monkey to rest on (**Figure 7B**) while they perform the
111 behavioral tasks.
112 The animals have access (**Figure 7B**) to a touch screen (e.g.,
113 android tablets, Apple iPads) via two apertures (**Figure 7C**) in the
114 front of the kiosk. The top aperture (see **Figure 7C** for dimensions)
115 provides a clear, non-occluded visual access to the screen. It
116 currently has a width of 6 inches which is broad enough to fit the
117 monkeys where the maximum distance between the edges of the
118 two eyes are usually \sim 3 inches. The bottom aperture allows
119 the monkey to reach inside with their hands and use the touch
120 screen. Of note, it is important to make sure that the edges of this
121 aperture have been additionally smoothed by the use of sand
122 paper (or something similar), since rubbing along those edges
123 (if not polished) might lead to minor cuts in the NHP arms. It
124 is usually advised that the experimenter tests the kiosk in a lab
125 setting with themselves (as subjects) first to ensure that everything
126 is working correctly. In the beginning of the NHP training pipeline,
127 the space between the touchscreen and the front screen may be
128 utilized to place toys and treats (**Figure 7D**) so that the monkey
129 feels motivated (and not afraid) to insert their hands through the
130 aperture.

131 **Concluding Remarks**

132 The kiosk presented here should be considered as a work in
133 progress. As we develop the system further, this article will be up-
134 dated with that information. We encourage the NHP neuroscience
135 research community to contact us with suggestions or clarification
136 questions. This will help us improve the system and make it more
137 versatile. Our next goals are to add a low cost eye tracking sys-
138 tem (as mentioned in **Figure 3B**), and to release the code-base
139 and material for the electronics that we are using. Our kiosks are
140 fully compatible with behavioral testing of macaques paired with
141 chemogenetic perturbations (currently used in this study: [Kar et](#)
142 [al., 2021](#)), and wireless electrophysiological recordings. As we
143 further develop those studies and techniques at York University,
144 we shall update the system details through this article (or other
145 similar documents).

146 **Open-source material**

147 The details of the kiosk can be found at the following GitHub page:
148 <https://github.com/vital-kolab/nhp-turk>

149 **Acknowledgments**

150 We thank Xiaogang Yan for assistance with cage compatibility
151 assessments. This work has been motivated by earlier work by
152 Elias Issa and other members of the DiCarlo Lab at the McGovern
153 Institute for Brain Research at MIT. We especially thank Sarah
154 Goulding for her help. The work has been supported by funds from
155 Canada Foundation for Innovation (CFI) and Simons Foundation
156 Autism Research Initiative (SFARI). This research was undertaken
157 thanks in part to funding from the Canada First Research Excel-
158 lence Fund (CFREF) administered through the VISTA Program at
159 York University. KK was supported by the Canada Research Chair
160 Program.

References

- | | |
|---|--|
| Rajalingham, R., Issa, E. B., Bashivan, P., Kar, K., Schmidt, K., & DiCarlo, J. J. (2018). Large-scale, high-resolution comparison of the core visual object recognition behavior of humans, monkeys, and state-of-the-art deep artificial neural networks. <i>Journal of Neuroscience</i> , 38(33), 7255-7269. | 161
162
163
164
165
166 |
| Kar, K., Kubilius, J., Schmidt, K., Issa, E. B., & DiCarlo, J. J. (2019). Evidence that recurrent circuits are critical to the ventral stream's execution of core object recognition behavior. <i>Nature neuroscience</i> , 22(6), 974-983. | 167
168
169
170 |
| Kar, K., Schrimpf, M., Schmidt, K., DiCarlo, J. J. (2021). Chemogenetic suppression of macaque V4 neurons produces retinotopically specific deficits in downstream IT neural activity patterns and core object recognition behavior. <i>Journal of Vision</i> , 21(9), 2489-2489. | 171
172
173
174 |